





GeneStream
DNA-based research tools


RapidReporter[®]

Example Data




Induction of CRE: Flasks of 293 cells were transiently transfected with the indicated plasmids, and 6h later were split evenly into 96-well plates such that all samples for a given plasmid were derived from the same transfection. After an overnight incubation, samples were either left untreated or received either 10 μ M forskolin or 4 μ M isoproterenol. Samples were then lysed at the indicated time-points and measured for luminescence using either the Rapid-G1 kit (for Ga plasmids), or a standard firefly kit (for FF plasmids). Data are plotted as mean +/- SEM.





Induction of NFκB: Flasks of 293 cells were transiently transfected with the indicated plasmids, and 6h later were split evenly into 96-well plates such that all samples for a given plasmid were derived from the same transfection. After an overnight incubation, samples were either left untreated or received 10ng/ml TNFα. Samples were then lysed at the indicated time-points and measured for luminescence using either the Rapid-G1 kit (for Ga plasmids), or a standard firefly kit (for FF plasmids). Data are plotted as mean +/- SEM.

Repression of TRE



Repression of TRE: Flasks of HeLa Tet-off cells were transiently transfected with the indicated plasmids, and 6h later were split evenly into 96-well plates such that all samples for a given plasmid were derived from the same transfection. After an overnight incubation, samples were either left untreated or received 3 μ g/ml doxycycline. Samples were then lysed at the indicated time-points and measured for luminescence using either the Rapid-G1 kit (for Ga plasmids), or a standard firefly kit (for FF plasmids). Data are plotted as mean \pm SEM.